



Effervescent Disinfectant Tablets



MAINTAIN A SAFE SPACE WITH ONE TABLET

Now, more than ever, it is crucial to keep a clean, safe and healthy environment.

Your commitment to maintaining a safe space can save lives.



SURIV, a hospital-grade, multi-purpose disinfectant and sanitizing agent. The active ingredient of **SURIV** is **NaDCC**, a disinfectant approved by the **US-EPA*** for **Coronavirus** and recommended by the **CDC***, **WHO*** and **UNICEF***.



Effervescent Disinfectant Tablets
200 NaDCC Tablets
Each tablet weighs 3.25g



Effervescent Disinfectant Extra Strong Tablets
200 NaDCC Tablets
Each tablet weighs 5g



Detergent Sanitizer Tablets
Cleans & Disinfects
40 NaDCC Tablets
Each tablet weighs 3.3g



Safe for users, pets, plants, clothing & surfaces



User-friendly, easy to use & exact



Resistant to spoilage with a long shelf life



Inexpensive & cost effective



Virucidal with a vicious killing power

3 EASY STEPS



Fill your container with water



Drop a tablet in your container



Use after tablet dissolves

For optimum virus disinfection a minimum of 1,000 ppm is required

RECOMMENDED DILUTION:

*Add 1.67g tablet to 1 liter water

*Add 2.5g tablet to 1.5 liters water

PUBLIC HEALTH SAFETY AND INFECTION CONTROL

Household surface sanitizer



Outbreak & biohazard disinfectant



General items sanitizer



Hospital surface & equipment disinfectant



Food service sanitizer




Public space disinfectant



Deodorizer



Laundry disinfectant



CHLORINE DISINFECTANT TABLETS

Suriv is especially formulated as a safe and easy to prepare solution with broad spectrum activity against microorganisms including Bacteria, Fungi, Spores and Viruses. **Suriv** is suitable for general sanitizing, deodorizing and disinfecting of most non-porous hard surfaces. **Directions for use:** Fill container with water. To determine how much water to fill, please refer to dilution table below for guidance. Drop one tablet into the container and wait for it to dissolve before use. Warm water will result in faster dilution.

SURIV 1.67g DILUTION TABLE

1 TABLET IN	GIVES AVAILABLE CHLORINATION OF	APPLICATIONS
100 mL	10,000 ppm	Body fluid spills & conditions of heavy soiling
400 mL	2,500 ppm	Non-metallic medical & veterinary appliances & laboratory discard jars
1 liter	1,000 ppm	General disinfection
2.5 liters	400 ppm	Toilets, drains & sinks
5 liters	200 ppm	Food preparation areas, non-metallic equipment, surfaces & materials
8 liters	125 ppm	Baby bottle sterilization
10 liters	100 ppm	Cloths & mops

SURIV XS 2.5g DILUTION TABLE

1 TABLET IN	GIVES AVAILABLE CHLORINATION OF	APPLICATIONS
150 mL	10,000 ppm	Body fluid spills & conditions of heavy soiling
600 mL	2,500 ppm	Non-metallic medical & veterinary appliances & laboratory discard jars
1.5 liters	1,000 ppm	General disinfection
3.75 liters	400 ppm	Toilets, drains & sinks
7.5 liters	200 ppm	Food preparation areas, non-metallic equipment, surfaces & materials
12 liters	125 ppm	Baby bottle sterilization
15 liters	100 ppm	Cloths & mops

SURIV DST 1.67g DILUTION TABLE

1 TABLET IN	GIVES AVAILABLE CHLORINATION OF	APPLICATIONS
100 mL	10,000 ppm	Body fluid spills & conditions of heavy soiling
400 mL	2,500 ppm	Non-metallic medical & veterinary appliances & laboratory discard jars
1 liter	1,000 ppm	General disinfection
3 tablets in 5 liters	600 ppm	Stainless steel instruments & materials
2.5 liters	400 ppm	Toilets, drains & sinks
5 liters	200 ppm	Food preparation areas, non-metallic equipment, surfaces & materials
10 liters	100 ppm	Cloths & mops



Do not take internally.



Keep out of reach of children.



Keep the container tightly closed in a cool dry place.



Do not inhale.



Avoid contact with eyes or skin.



Avoid prolonged contact with stainless steel equipments.



Do not mix with acid or cationic detergents